

CONVENTIONAL BARCODING

ONTbarcoder (available from <https://github.com/asrivathsan/ONTbarcoder>) is installed by unzipping the folder with the version of the program that supports the operating system on your computer. The folders are available from releases in Github. For MacOS, an .app bundle has been created.

ONTbarcoder is designed to handle dual tagged amplicon pools. It performs demultiplexing and consensus calling. If the software is used for the whole process, the format specifications are relatively simple and you can skip to the next section. If however, the demultiplexing is carried out independently, please follow the guidelines in Format Specifications, Point 2 in Page 10 of this documentation. It would be important to ensure that consensus calling is done under the assumption that primers and tags are removed. At least, tags should be removed if primers do not have ambiguous bases. The length parameters should be modified to expected consensus length.

ONTbarcoder has been optimized for **protein coding gene like COI**. While there are ways to obtain consensus for length variable non coding genes, this has not been extensively tested. *Please refer to last page of the manual to see which outputs to use if length variable and noncoding genes are being barcoded.*


ONTbarcoder's three modules.

1. Demultiplexing: MinION reads are assigned to specimen-specific bins.
2. Barcode calling: The reads in the specimen-specific bins are used to derive the barcodes based on alignment and consensus calling.
3. Barcode comparison module: Two or more sets of barcodes can be compared.


We recommend that you try the pipeline using our DatasetA containing 257 amplicons generated from Flongle. It is available from [here](#). Information to download the various datasets is available [here](#). Prior to running any dataset, ensure you have at least 1.5 times the space required for the input FASTQ file.

A [simple video tutorial](#) is available. This tutorial is for ONTbarcoder1 which describes Conventional barcoding. You can use "Conventional barcoding" of ONTbarcoder2, you don't need to download the older ONTbarcoder.

1. Open ONTbarcoder, select Conventional barcoding




This should lead to the following menu. The rest of the steps are identical to ONTbarcoder1 shown below. Only modification in ONTbarcoder2 is the option for selecting tag errors



Demultiplexing and barcode calling using the “Demultiplexing-Consensus calling tab”

2. Drag in Input Fastq and [demultiplexing file](#), either simultaneously or one by one




3. Press "Let's go"

4. Configure your run based on the settings described below. DatasetA does not require any changes to the default. The default assumes that you are working on invertebrate barcodes that are 658 bp in length. The remaining settings define which reads are demultiplexed, where the primers are found, how many reads are used in the different steps, and how the consensus barcodes are determined.

BASIC SETTINGS

Analysis workflow



Demultiplexing Settings

- Minimum Length: 658 (Annotation: Reads longer than this are accepted for demultiplexing)
- Length of barcode: 658 (Annotation: Expected barcode length, used for QC and determining expected product length)
- Window to define product length (length of barcode +/-): 100
- Window for primer and tag search: 100

Consensus by length settings

- Coverage used: 25,50,100,200,500 (Annotation: Number of reads used for consensus by length. If multiple coverages are specified (delimiter=comma), then iterative mode is used. In iterative mode, consensus are called at the lowest coverage. Those barcodes not passing the QC criteria of translation, meeting barcode length and being ambiguity free are passed for barcode calling at higher coverage.)
- Maximum deviation of read length from barcode length: 50

Consensus by similarity settings

- Coverage used: 100 (Annotation: Number of reads used for consensus by similarity.)

General settings for consensus calling

- Main consensus calling frequency: 0.3
- Range of frequencies to assess: 0.2,0.5
- to be examined at step size of: 0.05
- Genetic Code: 5. The Invertebrate Mitochondrial Code (Annotation: Genetic code)

Steps to run: ☒ Consensus by length ☒ Consensus by similarity ☒ Consensus by barcode comparisons

Run (Annotation: Start the run)

ADVANCED SETTINGS

The screenshot shows a software interface for 'ADVANCED SETTINGS' with several sections and input fields. Green arrows point from callout boxes to specific settings:


- Demultiplexing Settings**
 - Minimum Length:** 658
 - Length of barcode:** 658
 - Window to define product length (length of barcode +/-):** 100 (Callout: Window size for determining product length (explen+primer/tag length)+- window size determine product length)
 - Window for primer and tag search:** 100
- Consensus by length settings**
 - Coverage used:** 25,50,100,200,500
 - Maximum deviation of read length from barcode length:** 50
- Consensus by similarity settings**
 - Coverage used:** 100
- General settings for consensus calling**
 - Main consensus calling frequency:** 0.3 (Callout: Main consensus calling frequency, i.e. at least 30% of the bases should be occupied by a single nucleotide for a consensus to be called)
 - Range of frequencies to assess:** 0.2,0.5
 - to be examined at step size of:** 0.05
- Genetic Code:** 5. The Invertebrate Mitochondrial Code
- Steps to run:**
 - ☒ Consensus by length (Callout: Enables/disables **Consensus by length**: Step that uses reads closest to the known length of the barcode are used for consensus calling.)
 - ☒ Consensus by similarity (Callout: Enables/disables **Consensus by similarity**: Step that improves consensus in previous step by taking reads most similar to it.)
 - ☒ Consensus by barcode comparisons (Callout: Enables/disables **Consensus by barcode comparisons**, i.e. barcode fixing that fills/eliminates gaps in barcodes by comparison to other barcodes in the run)

Run


Callout boxes:

- Primers & tags are searched within the first/last specified number of bp of a read (points to 'Window for primer and tag search')
- Other consensus calling frequencies to test (points to 'Range of frequencies to assess')
- Specifies the steps for the range. In current settings 0.2,0.25,0.3,0.35,0.4,0.45,0.5 are tested (points to 'to be examined at step size of')


5. Select an output folder. If you select “default”, it will create a folder within the directory containing the software with a unique name starting with ONTbarcoderOutput. If you choose “Select Directory” ensure that the directory name has no empty space and is empty. Generally avoid complex characters in your names (spaces, brackets, slashes, colons)



6. The analysis starts, a progress bar appears, and results are generated in respective folders as the analysis advances.



7. The analysis is complete and all output files have been created.



8. Results Table: The table reports that 242 "Filtered" barcodes were obtained of which 192 are QC-Compliant. This means they have no ambiguities, are translatable, have expected barcode length and are not having any gaps in an internal Multiple Sequence Alignment (MSA) check. The filtered sets include addition 50 barcodes that have <1% ambiguous bases, are translatable, have expected length, and may have been corrected for up to 5 indels. The table can be copied but it also stored in the "[Outputtable.csv](#)" file in the output folder.


ONTBarcoder

Demultiplexing - Consensus calling Compare barcode sets Results of Barcode calling


	1	2
1 Number of reads in file		294896
2 Number of reads passing length filter		222189
3 Number of reads used for demultiplexing		190952
4 Number of samples in demultiplexing file		264
5 Number of samples with >=5X coverage		250
6 Number of good barcodes obtained after first alignment of up to 200 reads		167
7 Number of erroneous barcodes obtained after first alignment of up to 200 r...		69
8 Number of good barcodes obtained after aligning similar reads		15
9 Number of erroneous barcodes obtained after aligning similar reads		68
10 Number of barcodes fixed		69
11 Final number of barcodes of expected length		251
12 Final number of barcodes that cannot be fixed		0
13 Number of Ns in final barcodes		176
14 Number of filtered barcodes		242
15 Number of QC_Compliant barcodes		192
16 Number of barcodes with 1-5 errors		48
17 Number of barcodes with 6-10 errors		5
18 Number of barcodes with 11-15 errors		1

Quit

9. The output folder contains the following



10. The Main Barcode Results folder contains the main barcode output file "**Filtered_barcodes**" fasta file. The barcodes without ambiguous bases are in the QC_Compliant_barcodes file. The zip folders contain demultiplexed reads per dataset divided into each of the categories. All barcodes should be checked via BLAST for contamination and/or verification that the barcode obtained belongs to the expected taxon.



Summary

Expected input and output

1. Conventional barcoding	Input required	Results
Mode 1: Demultiplexing + barcode calling	<ul style="list-style-type: none"> FASTQ file obtained after base-calling Demultiplexing file 	<ul style="list-style-type: none"> Demultiplexed reads in “demultiplexed” folder Overall summary in “Outputtable.csv” Per barcode summary in “runsummary.xlsx”
Mode 2: Barcode calling only: Barcodes are derived from reads in specimen-specific bins	<ul style="list-style-type: none"> Folder containing Demultiplexed FASTA files. Format specified under “Directory containing Demultiplexed FASTA files” 	<ul style="list-style-type: none"> Barcode sets in “Main barcode results” folder Overall summary in “Outputtable.csv” Per barcode summary in “runsummary.xlsx”
Mode 3: Improving barcodes using “Consensus by similarity” and “Consensus by barcode comparisons”	<ul style="list-style-type: none"> Sequences in “barcodesets” folder 	<ul style="list-style-type: none"> Barcode sets available from “Main barcode results” folder Overall summary in “Outputtable.csv” Per barcode summary in “runsummary.xlsx”

Format Specifications

File	
1. FASTQ file	Standard fastq, generated after basecalling with ONT software
2. Demultiplexing File	<p>A 5-column csv file with the following headers: SpecimenID, TagFsequence, TagRsequence, PrimerF, PrimerR</p> <p>You can only demultiplex one one primer pair at a time. FASTQ files with data for multiple pairs, have to processed sequentially.</p> <p>Please avoid unusual characters in Specimen ID (e.g. characters like “(){}[]V.,;*\$” will lead to crashes)</p>
3. Directory containing Demultiplexed FASTA files	Input for Mode 2 should be a directory that contains only FASTA files. Directory should not be empty nor should it contain other files. The fasta files must have names as “sampleID_all.fa” i.e. the suffix _all.fa is critical
4. Input file for improvement (Step 4)	ONTbarcoder’s pipeline can be started at different points. Demultiplexing files can be used to only carry out barcode calling. Consensus by barcode can be started with files that have the format specified in the hyperlink

LENGTH VARIABLE AND NON-CODING GENES

Accurate barcodes have been obtained by preselection of reads to the specified length criteria of the gene. Barring this issue, ONTbarcoder can be used for non-coding, length variable genes. For this, please disable consensus by similarity and consensus by barcode fixing step in the menu. The barcode calling can be conducted as per normal. The resulting barcode file can be found in barcodesets/consensus_all_step1.fa

This can be parsed by the user to retain barcodes with few ambiguities (for e.g. <1%)